

FibroScan[®]
compact 530

ROLL-STAND



PROVIDING A NEW
LEVEL OF SUPPORT



FibroScan[®] compact 530 *ROLL-STAND*

As you move your FibroScan 530 Compact among locations, the Roll-stand can provide a sturdy and functional platform that further protects your investment.

COMFORT AND CONVENIENCE

- Adjust the height with the touch of a lever
- Probe and gel holders are easily removed for cleaning
- Shelf for adding an optional keyboard

SECURITY

- Lockable drawer to store and secure the probes


MOBILITY

- Moves easily with four locking wheels

FLEXIBILITY

- Secure installation is quick and easy (only 3 screws)

RAISED HEIGHT POSITION:	43.5"
LOWER HEIGHT POSITION:	34"
PRINTER SHELF DIMENSIONS:	15.5" x 12"
PRINTER SHELF MAXIMUM LOAD:	22 lbs.
WEIGHT (EMPTY):	83.6 lbs.



Visit us at www.echosens.us

ECHOSSENS NORTH AMERICA • 1050 Winter St. • Waltham, MA 02451
(781) 790-0845 • inquiries@echosens.com